

NEW Agilent MSD Vacuum Pump (GLD-040)



NEW Vacuum Pump for Agilent MSD (GLD-040)

Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [Agilent](#)

Description

NEW Vacuum Pump for Agilent MSD (GLD-040)

The GLD series rotary vane pumps from ULVAC combine high performance, low vibration, and low noise. Design features include: 1.7 cfm; gas ballasting to remove solvent vapors; forced oil supply for optimum lubrication during pump-down; and an anti-suckback device to protect the chamber during power failure.

Includes:

- **NEW ULVAC Vacuum Pump (GLD-040)**
- **ONE (1) YEAR WARRANTY.**

Specifications:

Stages	2
Pumping Speed (cfm)	1.7

Ultimate pressure (Torr) (using McLeod gauge)	5x10 ⁻⁴
Voltage (50/60 Hz)	100-115/200-230 V
Inlet	KF25
Exhaust	KF25
Weight	14 kg (31lbs)
Dimensions (mm) (W x L x H)	165 x 399 x 223